

## Anti-DIO2 antibody [R06-4G-8] (STJA0033128)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	WB 1:500-1:1000 IHC-P 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG2ak
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1734</a>
Gene Symbol	<a href="#">DIO2</a>
UniProt ID	<a href="#">IOD2_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---