

Anti-OMA1 antibody [R05-8B-2] (STJA0033131)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC-P 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG1k
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	115209
Gene Symbol	OMA1
UniProt ID	OMA1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



OMA1 antibody immunohistochemistry mouse liver STJA0033131