

**Anti-BRD8 antibody [R09-7X-7] (STJA0033137)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P           |
| Host / Source | Mouse              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                           |
|---------------------|-------------------------------------------|
| Clonality           | Monoclonal                                |
| Conjugation         | Unconjugated                              |
| Purification        | Protein G                                 |
| Dilution Range      | WB 1:200-1:1000<br>IHC-P 1:50-1:100       |
| Formulation         | Liquid in PBS, Glycerol and BSA           |
| Isotype             | IgG2ak                                    |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">10902</a>      |
| Gene Symbol | <a href="#">BRD8</a>       |
| UniProt ID  | <a href="#">BRD8_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|