

Anti-EIF3L antibody [R04-5C-5] (STJA0033141)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	WB 1:2000-1:10000
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG2Ak
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	51386
Gene Symbol	EIF3L
UniProt ID	EIF3L_HUMAN

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



EIF3L antibody western blot 293 Raw2647 3T3 HepG2 STJA0033141