

Anti-CD5 antibody [R05-4K-9] (STJA0033145)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IF/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	ICC:IF 1:50 FC 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG2ak
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	921
Gene Symbol	CD5
UniProt ID	CD5_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------