

Anti-CD33 antibody [R07-3M-6] (STJA0033147)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	IHC-P 1:200-1:1000
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG2ak
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	945
Gene Symbol	CD33
UniProt ID	CD33_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CD33 antibody western blot Raji Daudi Ramos THP-1 STJA0033147