

Anti-CD3 epsilon antibody (STJA0033158)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A affinity column
Dilution Range	IHC-F-1:50-1:100 FC-1:50-1:100
Formulation	Liquid in PBS, pH 7.4
Isotype	IgG2a
Storage Instruction	Store at 2-8°C. Do not freeze.

TARGET INFORMATION

Gene ID	916
Gene Symbol	CD3E
UniProt ID	CD3E_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---