

Anti-alpha 1 Fetoprotein antibody (STJA0033159)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/ELISA
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	affinity purified by Protein G
Dilution Range	WB-1:500-1:2000 IHC-1:100-1:500 IF-1:100-1:500 ELISA-1:5000-1:10000
Formulation	Liquid in 0.01M PBS, PH7.4, with 0.02% Proclin300.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	174
Gene Symbol	AFP
UniProt ID	FETA_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---