

## Anti-eNOS antibody (STJA0033161)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	affinity purified by Protein G
Dilution Range	WB-1:500-1:1000 IHC-1:100-1:500 IF-1:100-1:500
Formulation	Liquid in 0.01M TBS, PH7.4, with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">4846</a>
Gene Symbol	<a href="#">NOS3</a>
UniProt ID	<a href="#">NOS3_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---