

Anti-SPHK2 antibody (STJA0033165)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/ICC/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	affinity purified by Protein A
Dilution Range	IHC-1:400-1:800 IF-1:100-1:500
Formulation	Liquid in 0.01M TBS, PH7.4, with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	56848
Gene Symbol	SPHK2
UniProt ID	SPHK2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---