

Anti-beta Actin antibody [8F10-G10] {HRP} (STJA0033171)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Mouse/Monkey/Goat/Rat/Hamster/S.cerevisiae/Arabidopsis/P.pastoris/Rabbit/Zebrafish/D.melanogaster/Rice Leaf/Cucumber/C.reinhardtii

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2.4 mg/mL
Conjugation	HRP
Purification	Affinity Purified
Dilution Range	WB 1:5000-1:10000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	60
Gene Symbol	ACTB
UniProt ID	ACTB_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



beta actin antibody western blot HeLa 3T3 C6 CHO-K1 COS7 STJA0033171



Beta-actin antibody western blot rice leaf STJA0033171