

Anti-beta Actin antibody [6D7-F6-G10] (STJA0033177)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Rat/Mouse/Goat/Monkey/Hamster

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2.61 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:1000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG3
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	60
Gene Symbol	ACTB
UniProt ID	ACTB_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Beta-actin antibody western blot HeLa 3T3 COS7
 CHO-K1 goat muscle STJA0033177