

Anti-GAPDH antibody [7E4-H6-H6] (STJA0033182)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Mouse
Reactivity	Human/Monkey/Rat/Mouse/S.cerevisiae/P.pastoris/Hamster

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1.9 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:5000-1:10000
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgM
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	2597
Gene Symbol	GAPDH
UniProt ID	G3P_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



GAPDH antibody western blot HeLa 3T3 C6 COS7 CHO-K1 STJA0033182



GAPDH antibody western blot S. cerevisiae P. pastoris STJA0033182