

Anti-DNA PKCs antibody [8D3-B12-F11] (STJA0033202)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5591
Gene Symbol	PRKDC
UniProt ID	PRKDC_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




DNAPKcs antibody western blot HeLa K562 HL-60 MOLT4 STJA0033202



DNA-PKcs antibody immunocytochemistry HeLa cells 8D3 STJA0033202



DNA-PKcs antibody immunoprecipitation HeLa lysates STJA0033202



DNA-PKcs antibody western blot K562 Molt4 A549 HCT116 8D3 STJA0033202