

Anti-EGFR antibody [6H11-2D11-G3] (STJA0033225)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF/IP
Host / Source	Mouse
Reactivity	Human/Monkey

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.45 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	1956
Gene Symbol	EGFR
UniProt ID	EGFR_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




EGFR antibody western blot A549 COS7 lysates STJA0033225



EGFR antibody immunocytochemistry HeLa cells 6H11 STJA0033225



EGFR antibody immunoprecipitation HeLa lysates 6H11 STJA0033225



EGFR antibody immunocytochemistry MDA-MB-468 cells STJA0033225



EGFR antibody western blot A549 MDA-MB-468 lysates STJA0033225