

Anti-Flag-Tag antibody [3A6-E5-C10] (STJA0033233)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Mouse
Reactivity	Species-independent

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	3.6 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:1000-1:10000 IP 1:20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



FLAG Tag antibody western blot CHO-K1 cells
STJA0033233



Flag Tag antibody immunoprecipitation CHO K1
2011263A6 STJA0033233