

**Anti-Clostridium Difficile Toxin A antibody [2D4-D3] (STJA0033274)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA              |
| Host / Source | Mouse              |
| Reactivity    | C.difficile        |

**PRODUCT PROPERTIES**

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | ELISA 1:10000                                                                            |
| Formulation         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.          |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|