

## Anti-Thymidine Phosphorylase antibody [ZN98] (STJA0033311)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.5-1.0 mg/mL                                                                            |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Chromatography                                                                  |
| Dilution Range      | IHC 1:50-1:100                                                                           |
| Formulation         | PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                     |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">1890</a>       |
| Gene Symbol | <a href="#">TYMP</a>       |
| UniProt ID  | <a href="#">TYPH_HUMAN</a> |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Thymidine Phosphorylase antibody  
immunohistochemistry human lymph node  
STJA0033311



Thymidine Phosphorylase antibody  
immunohistochemistry human colon  
adenocarcinoma STJA0033311