

## Anti-SLC7A8 antibody [ZN70] (STJA0033336)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.5-1.0 mg/mL                                                                            |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Chromatography                                                                  |
| Dilution Range      | IHC 1:50-1:100                                                                           |
| Formulation         | PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                     |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">23428</a>      |
| Gene Symbol | <a href="#">SLC7A8</a>     |
| UniProt ID  | <a href="#">LAT2_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SLC7A8 antibody immunohistochemistry human colon tissue STJA0033336



SLC7A8 antibody immunohistochemistry human kidney tissue STJA0033336