

## Anti-PBX1 antibody [ZN110] (STJA0033355)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.5-1.0 mg/mL                                                                            |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Chromatography                                                                  |
| Dilution Range      | IHC 1:50-1:100                                                                           |
| Formulation         | PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                     |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">5087</a>       |
| Gene Symbol | <a href="#">PBX1</a>       |
| UniProt ID  | <a href="#">PBX1_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



PBX1 antibody immunohistochemistry human lung carcinoma STJA0033355



PBX1 antibody immunohistochemistry human bladder carcinoma STJA0033355