

## Anti-IDO1 antibody [ZN126] (STJA0033368)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.5-1.0 mg/mL                                                                            |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Chromatography                                                                  |
| Dilution Range      | IHC 1:50-1:100                                                                           |
| Formulation         | PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                     |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">3620</a>        |
| Gene Symbol | <a href="#">IDO1</a>        |
| UniProt ID  | <a href="#">I23O1_HUMAN</a> |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IDO1 antibody immunohistochemistry human lung carcinoma STJA0033368



IDO1 antibody immunohistochemistry human placenta STJA0033368