

Anti-BOB1 antibody [4G9-4A8-3E8] (STJA0033414)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	IHC 1:100-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2a
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5450
Gene Symbol	POU2AF1
UniProt ID	QBF1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



BOB1 antibody immunohistochemistry Hodgkin lymphoma STJA0033414



BOB1 antibody immunohistochemistry Hodgkin lymphoma paraffin-embedded STJA0033414



BOB1 antibody immunohistochemistry paraffin-embedded tonsil tissue STJA0033414