

Anti-Calponin 1 antibody [9E8-1H9-8C6] (STJA0033429)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	IHC 1:100-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	1264
Gene Symbol	CNN1
UniProt ID	CNN1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Calponin 1 antibody immunohistochemistry laryngeal squamous cell carcinoma STJA0033429



Calponin 1 antibody immunohistochemistry paraffin-embedded leiomyoma STJA0033429



Calponin 1 antibody immunohistochemistry cervical squamous carcinoma STJA0033429