

## Anti-PSMA antibody [5D5-7G4-1B7] (STJA0033446)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 1 mg/mL                                                                                  |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Chromatography                                                                  |
| Dilution Range      | IHC 1:100-1:500                                                                          |
| Formulation         | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.          |
| Isotype             | IgG1                                                                                     |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">2346</a>        |
| Gene Symbol | <a href="#">FOLH1</a>       |
| UniProt ID  | <a href="#">FOLH1_HUMAN</a> |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



PSMA antibody immunohistochemistry prostate tissue STJA0033446



PSMA antibody immunohistochemistry paraffin-embedded prostate tissue STJA0033446



PSMA antibody immunohistochemistry prostate tissue paraffin-embedded STJA0033446