

Anti-p57 Kip2 antibody [7H9-2G7-5G5] (STJA0033473)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	IHC 1:100-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG3
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	1028
Gene Symbol	CDKN1C
UniProt ID	CDN1C_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



p57 Kip2 antibody immunohistochemistry partial hydatidiform mole STJA0033473



p57 Kip2 antibody immunohistochemistry placenta paraffin-embedded STJA0033473



p57 Kip2 antibody immunohistochemistry placenta paraffin-embedded STJA0033473