

Anti-CD68 antibody [1B8-8F6-6H8] (STJA0033483)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.38mg/mL
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	IHC 1:100-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	968
Gene Symbol	CD68
UniProt ID	CD68_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




CD68 antibody immunohistochemistry acute myeloid leukemia paraffin-embedded STJA0033483



CD68 antibody immunohistochemistry lymphoma paraffin-embedded tissue STJA0033483



CD68 antibody immunohistochemistry malignant fibrous histiocytoma STJA0033483



CD68 antibody immunohistochemistry tonsil tissue paraffin-embedded STJA0033483