

Anti-Actin pan antibody [8H5-1H4-9E8] (STJA0033504)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	IHC 1:100-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG2b
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID
[71](#)
[60](#)
[NA](#)
[ACTG1](#)
[ACTB](#)
[ACTG_HUMAN](#)
[ACTB_HUMAN](#)
[ACTBM_HUMAN](#)

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Actin pan antibody immunohistochemistry cardiac muscle paraffin-embedded STJA0033504



Actin pan antibody immunohistochemistry colon tissue paraffin-embedded STJA0033504



Actin pan antibody immunohistochemistry skeletal muscle paraffin-embedded STJA0033504