

Anti-Collagen III alpha 1 antibody [6H3-9G8-1C4] (STJA0033514)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	IHC 1:100-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Isotype	IgG1
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	1281
Gene Symbol	COL3A1
UniProt ID	CO3A1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Collagen III alpha 1 antibody immunohistochemistry skin paraffin-embedded STJA0033514



Collagen III alpha 1 antibody immunohistochemistry skin paraffin-embedded STJA0033514



Collagen III alpha 1 antibody immunohistochemistry skin paraffin-embedded STJA0033514