

## Anti-ALK-R1275Q Mutant antibody (STJA0033556)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P/ELISA     |
| Host / Source | Mouse              |
| Reactivity    | Vertebrates        |

### PRODUCT PROPERTIES

|                     |                                                                           |
|---------------------|---------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                |
| Conjugation         | Unconjugated                                                              |
| Purification        | Purified from ascites                                                     |
| Dilution Range      | ELISA-1:1000-1:5000<br>WB-1:100-1:1000<br>IHC-P-1:50-1:100                |
| Formulation         | Liquid in PBS (without Mg2+ and Ca2+) , pH 7.4, 150 mM NaCl, 50% glycerol |
| Isotype             | IgG                                                                       |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles.                                 |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



ALKR1275Q antibody western blot recombinant protein STJA0033556