

## Anti-BRaf-N581S Mutant antibody (STJA0033562)

### GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/ICC/IF/IP/ELISA
Host / Source	Mouse
Reactivity	Vertebrates

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified from ascites
Dilution Range	ELISA-1:1000-1:5000 WB-1:500-1:2000 IHC-P-1:100-200 IF-1:50-1:200
Formulation	Liquid in PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ) , pH 7.4, 150 mM NaCl, 50% glycerol
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



BRAF N581S antibody western blot recombinant proteins STJA0033562



B-Raf antibody immunofluorescence HEK293T cells STJA0033562