

## Anti-CTNNB1-S45Y Mutant antibody (STJA0033581)

### GENERAL INFORMATION


|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ICC/IF/ELISA    |
| Host / Source | Mouse              |
| Reactivity    | Vertebrates        |

### PRODUCT PROPERTIES


|                     |                                                                             |
|---------------------|-----------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                  |
| Conjugation         | Unconjugated                                                                |
| Purification        | Purified from ascites                                                       |
| Dilution Range      | WB-1:500-1:1000<br>IHC-P-1:50-1:100<br>ELISA-1:1000-1:2000<br>IF-1:50-1:200 |
| Formulation         | Liquid in PBS (without Mg2+ and Ca2+) , pH 7.4, 150 mM NaCl, 50% glycerol   |
| Isotype             | IgG                                                                         |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles.                                   |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CTNNB1 S45Y antibody western blot recombinant protein STJA0033581



CTNNB1 antibody immunofluorescence HEK293T cells STJA0033581