

Anti-ErbB2-V777A Mutant antibody (STJA0033585)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/ELISA
Host / Source	Mouse
Reactivity	Vertebrates

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified from ascites
Dilution Range	WB-1:200-1:2000 ICC/IF-1:100 IHC-P-1:50-1:10 ELISA-1:500-1:5000
Formulation	Liquid in PBS (without Mg2+ and Ca2+) , pH 7.4, 150 mM NaCl, 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti-ErbB2V777A antibody immunofluorescence
 HEK293T cells STJA0033585