

Anti-EGFR-G719S Mutant antibody (STJA0033587)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/ELISA
Host / Source	Mouse
Reactivity	Vertebrates

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified from ascites
Dilution Range	ELISA-1:1000-1:5000 WB-1:1000-1:2000 IHC-P-1:50-1:1000 ICC/IF-1:50-1:100
Formulation	Liquid in PBS (without Mg ²⁺ and Ca ²⁺) , pH 7.4, 150 mM NaCl, 50% glycerol
Isotype	IgG2b
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

	
EGFRG719S antibody western blot recombinant protein STJA0033587	EGFR G719S antibody immunofluorescence HEK293T cells STJA0033587