

## Anti-EGFR-D761Y Mutant antibody (STJA0033592)

### GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/ELISA
Host / Source	Mouse
Reactivity	Vertebrates

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified from ascites
Dilution Range	WB-1:500-1:2000 IHC-P-1:50-1:100 ELISA-1:1000-1:5000 IF-1:50-1:200
Formulation	Liquid in PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ) , pH 7.4, 150 mM NaCl, 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



EGFRD761Y antibody western blot recombinant protein STJA0033592



EGFRD761Y antibody immunofluorescence HEK293T cells STJA0033592