

## Anti-EGFR-S768I Mutant antibody (STJA0033601)

### GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/ELISA
Host / Source	Mouse
Reactivity	Vertebrates

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified from ascites
Dilution Range	ELISA-1:1000-1:5000 WB-1:1000-1:5000 IHC-P-1:50-1:200 IF-1:50-1:200
Formulation	Liquid in PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ) , pH 7.4, 150 mM NaCl, 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



EGFR S768I antibody western blot recombinant protein STJA0033601



EGFRS768I antibody immunofluorescence HEK293T cells STJA0033601