

## Anti-IDH1-R132P Mutant antibody (STJA0033621)

### GENERAL INFORMATION


|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IHC-P/ICC/IF/ELISA |
| Host / Source | Mouse                 |
| Reactivity    | Vertebrates           |

### PRODUCT PROPERTIES


|                     |                                                                                 |
|---------------------|---------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                      |
| Conjugation         | Unconjugated                                                                    |
| Purification        | Purified from ascites                                                           |
| Dilution Range      | WB-1:100-1:1000<br>IHC-P-1:50-1:100<br>ICC/IF-1:50-1:100<br>ELISA-1:1000-1:5000 |
| Formulation         | Liquid in PBS (without Mg2+ and Ca2+) , pH 7.4, 150 mM NaCl, 50% glycerol       |
| Isotype             | IgG                                                                             |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles.                                       |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IDH1R132P antibody western blot recombinant proteins STJA0033621



IDH1R132P antibody immunofluorescence HEK293T cells STJA0033621