

## Anti-IDH2-R172S Mutant antibody (STJA0033622)

### GENERAL INFORMATION


|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IHC-P/ICC/IF/ELISA |
| Host / Source | Mouse                 |
| Reactivity    | Vertebrates           |

### PRODUCT PROPERTIES


|                     |                                                                               |
|---------------------|-------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                    |
| Conjugation         | Unconjugated                                                                  |
| Purification        | Purified from ascites                                                         |
| Dilution Range      | WB-1:100-1:1000<br>ICC/IF-1:50-1:100<br>ELISA-1:1000-1:5000<br>IHC-1:50-1:100 |
| Formulation         | Liquid in PBS (without Mg2+ and Ca2+) , pH 7.4, 150 mM NaCl, 50% glycerol     |
| Isotype             | IgG                                                                           |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles.                                     |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IDH2R172S antibody western blot recombinant protein STJA0033622



IDH2 antibody immunofluorescence HEK293T cells STJA0033622