

## Anti-KIT-W557R Mutant antibody (STJA0033627)

### GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/ELISA
Host / Source	Mouse
Reactivity	Vertebrates

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified from ascites
Dilution Range	ELISA-1:1000-1:5000 WB-1:500-1:1000 ICC/IF-1:50-1:100 IHC-P-1:25-1:100
Formulation	Liquid in PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ) , pH 7.4, 150 mM NaCl, 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



KITW557R antibody western blot recombinant KIT W557R antibody immunofluorescence protein STJA0033627 HEK293T cells STJA0033627



KIT W557R antibody immunofluorescence HEK293T cells STJA0033627