

Anti-KIT-V559D Mutant antibody (STJA0033628)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	WB/ICC/IF/ELISA
Host / Source	Mouse
Reactivity	Vertebrates

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified from ascites
Dilution Range	WB-1:500-1:1000 IHC-P-1:50-1:100 ELISA-1:1000-1:5000 IF-1:50-1:200
Formulation	Liquid in PBS (without Mg ²⁺ and Ca ²⁺) , pH 7.4, 150 mM NaCl, 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



KITV559D antibody immunofluorescence
HEK293T cells STJA0033628