

## Anti-Beta-Catenin-E53K Mutant antibody (STJA0033677)

### GENERAL INFORMATION


|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IHC-P/ICC/IF/ELISA |
| Host / Source | Mouse                 |
| Reactivity    | Vertebrates           |

### PRODUCT PROPERTIES


|                     |                                                                             |
|---------------------|-----------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                  |
| Conjugation         | Unconjugated                                                                |
| Purification        | Purified from ascites                                                       |
| Dilution Range      | WB-1:500-1:2000<br>IHC-P-1:50-1:100<br>ELISA-1:1000-1:5000<br>IF-1:50-1:200 |
| Formulation         | Liquid in PBS (without Mg2+ and Ca2+) , pH 7.4, 150 mM NaCl, 50% glycerol   |
| Isotype             | IgG                                                                         |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles.                                   |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Beta-cateninE53K antibody western blot recombinant protein STJA0033677



Beta-catenin antibody immunofluorescence HEK293T cells STJA0033677