

Anti-2-Hydroxyisobutyryl-Histone H3-Lys14 antibody [1P9-B4-U2] (STJA0033689)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P
Host / Source	Mouse
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC-P 1:50-1:100
Formulation	Mouse IgG in phosphate buffered saline, glycerol and BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
TARGET INFORMATION	
Gene ID	3017/8339/8343/8344/8346/8347
Gene Symbol	H2BC4.H2BC6.H2BC7.H2BC8.H2BC10
UniProt ID	H2B1C_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.