

Anti-Phospho-RIPK1-S166 antibody (STJA0033693)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/FC/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	ProA affinity purified IgG
Dilution Range	WB 1:1000-1:2000 IP 1:50-1:500 ICC:IF 1:100-1:500 IHC-P 1:100-1:500 IHC-F 1:100-1:500
Formulation	Liquid in PBS 59%, Sodium azide 0.01%, Glycerol 40%, BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	19766
Gene Symbol	Ripk1
UniProt ID	RIPK1_MOUSE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



RIPK1 antibody western blot immunoprecipitation
 661W cells STJA0033693