

Anti-CD68 antibody [R02-6T-1] (STJA0033698)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/FC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A purified
Dilution Range	IHC-P 1:200-1:1000 FC 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	968
Gene Symbol	CD68
UniProt ID	CD68_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CD68 antibody immunohistochemistry human
cervical carcinoma STJA0033698