

**Anti-C4.4A/LYPD3 antibody [R05-6K-7] (STJA0033716)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	WB 1:500-1:2000 IHC-P 1:200-1:1000 ICC:IF 1:50
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">27076</a>
Gene Symbol	<a href="#">LYPD3</a>
UniProt ID	<a href="#">LYPD3_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---