

Anti-Thymidylate Synthase antibody [R01-5S-3] (STJA0033717)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB |
| Host / Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES


| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein A |
| Dilution Range | WB 1:500-1:2000 |
| Formulation | Liquid in PBS, Glycerol and BSA |
| Isotype | IgG |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 7298 |
| Gene Symbol | TYMS |
| UniProt ID | TYSY_HUMAN |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Thymidylate Synthase antibody western blot HeLa
Jurkat HepG2 RAW 2647 C6 STJA0033717