

## Anti-CD42b antibody [R03-5M-2] (STJA0033719)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC-P 1:200-1:1000
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG
Storage Instruction	Store at-20°C. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2811</a>
Gene Symbol	<a href="#">GP1BA</a>
UniProt ID	<a href="#">GP1BA_HUMAN</a>

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



CD42b antibody immunohistochemistry mouse spleen STJA0033719