

Anti-CBP antibody [R06-7T-9] (STJA0033724)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	WB 1:1000-1:3000 ICC:IF 1:50
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgGk
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	1387
Gene Symbol	CREBBP
UniProt ID	CBP_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---