

Anti-CD16a antibody [R08-3J-8] (STJA0033734)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-P |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein A |
| Dilution Range | IHC-P 1:200-1:1000 |
| Formulation | Liquid in PBS, Glycerol and BSA |
| Isotype | IgG |
| Storage Instruction | Store at-20°C. Avoid freeze/thaw cycles. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 2214 |
| Gene Symbol | FCGR3A |
| UniProt ID | FCG3A_HUMAN |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CD16a antibody immunohistochemistry human liver STJA0033734