

## Anti-ALOX15 antibody [R02-9X-5] (STJA0033737)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/ICC/IF/FC
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC-P 1:50-1:100 ICC:IF 1:50 FC 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">246</a>
Gene Symbol	<a href="#">ALOX15</a>
UniProt ID	<a href="#">LOX15_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---