

Anti-CDK8 antibody [R03-4H-9] (STJA0033743)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-P |
| Host / Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES


| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein A |
| Dilution Range | IHC-P 1:50-1:100 |
| Formulation | Liquid in PBS, Glycerol and BSA |
| Isotype | IgG |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 1024 |
| Gene Symbol | CDK8 |
| UniProt ID | CDK8_HUMAN |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CDK8 antibody immunohistochemistry mouse colon STJA0033743